


## DISCUS DMC.Connect App

Import DICOM images, videos, and data via QR-Code

### Simple import of images and reports directly into PACS

Imagine being able to import patients' examination data into your PACS quickly and easily via a QR code, without needing to connect workstations to the internet or having to manually enter annoying login credentials. Additionally, the security of your institution remains intact, and your staff is relieved from extra work.



### The solution is called DISCUS DMC.Connect App.

The DISCUS DICOM Media Center with QR code scanner enables efficient and secure import of DICOM images and reports from web portals and local network directories. The connection to the internet is made exclusively through a hardened Linux server running on a virtual machine. This ensures that neither the employee nor the workstation has direct internet access. With this system, you increase security and efficiency when importing data into your network and PACS.

### Benefit from many other advantages:

- Import from all existing portals on the market
- Securely import patients' examination data into the local PACS without an internet connection at the workstations
- Access to portals only through a secure, hardened Linux server
- Secure and fast QR code integration using supplied scanners for all workstations for patient registration
- No time-consuming manual entry of login credentials
- Automated data import into PACS
- Minimize errors by eliminating manual input
- Verification of patient data for DICOM compliance and completeness
  - unauthorized files are ignored
- Simple and intuitive operation without extensive training
- Support for importing all common DICOM objects to ensure maximum compatibility

Call us or visit our new website if you have any further questions.

+49 06222-938830

[www.dicom-disc.de](http://www.dicom-disc.de)

[discus@adr-ag.de](mailto:discus@adr-ag.de)